


F560

maccura

AUTOMATIC HEMATOLOGY ANALYZER


- 3rd generation technology
- 60 samples per hour
- Reliable performance for aging blood/abnormal sample




Advanced Technology

3rd generation Tech

Fluorescence staining to Nucleic Acid



Special fluorescent staining solution will dye DNA or RNA blandly while 2nd Generation chemistry staining reagents will dye Enzymes/particles in cytoplasm. we know that different cell has different concentration of DNA or RNA , which cause the depth of dying is different. the more DNA or RNA , the stronger fluorescent signal. Since the nucleic acid is the most specific part of cell, so the 3rd Generation is more sensitive to distinguish different leukocyte, especially the abnormal cells




Combine 3rd Generation technology with flow cytometry, A single-cell stream quickly passes through a channel in the middle, and every passing cell is detected by three beams of light from three directions to get size, granularity and nucleic acid information

FSL (Forward Scattered Light) mainly reflects the size of the cells, SSC (Side Scattered Light) mainly reflects size and number of particle in cells SFL (Side Fluorescence Light) mainly reflects the concentration of nucleic acid

Excellent performance





High sensitive to abnormal cells



Atypical lymphocyte and immature granular cell have strong nucleic acid fluorescent signal, after fluorescent staining, they are easier to be detected

Help to distinguish abnormal myeloid and gonorrhoea cells

Trustable performance




good correlation with comparison system

Powerful Functions


Comprehensive flag information

- clinical flag
 - (1) enhanced abnormal cell detection capacity
 - (2) Help diagnosis such as hypochromia anemia, neutropenia, etc
- maintenance flag
 - (1) powerful debug function
 - (2) one click to remove error



High efficiency

- Through put 60samples/hour




Multi-channels

- Independent Baso test channel
- Specid DIFF channel with blood
- Individual RBC/PLT channel



Test options

- Mode :CBC、CBC+DIFF
- Sample type : whole blood, capillary blood, pre-dilution blood
- Auto sample dilution





Easy-to-use software

- Simple daily operation
 - (1)Visual and intuitive software interface
 - (2)Convenient data management
- Easy maintenance
 - (1)one click to remove clog
 - (2)powerful debug functions



Smart graphical analysis



3D scattergram for accurate WBC differentiation and pathological sample flag

Histograms for WBC/RBC/PLT

SPECIFICATIONS

Measurement principle	Nucleic Acid Fluorescence Staining and Flow Cytometry
Measurement Parameter	24 report parameters (WBC,RBC,HGB,MCV,MCH,MCH-C, RDW-CW,RDW-SD,HCT,PLT,MPV,P-DW,PCT,P-LCR,BASO#,BASO%,NEUT#,NEUT%E0#,E0%,LYMPH#,LYMPH%,MONO#,MONO%) 4 Research Parameter (IG#,IG%, OTHER#,OTHER%) 4 Graphs (2D and 3D scattergram analysis, 3 histograms)
Throughput	60 T/H
Test Mode	CBC / CBC+DIFF
Sample Type	Whole blood / Capillary blood / Pre-dilution blood
Sampling Method	Manual sampling
Sample Volume	20ul
Reagent	GD-5 (Diluent) LH-5 (HGB Lyse) LD-5 (DIFF Lyse) DD-5 (Dye) CC-5(Clean solution)
Power requirement	100-240V \leq 250VA, 50/60Hz
Dimensions	550x700x600mm
Weight	55kg

LINEARITY RANGE

Parameter	Linear Measurement Range	Linear Tolerance	r
WBC	$1.0 \times 10^9/L \sim 10.0 \times 10^9/L$ $10.1 \times 10^9/L \sim 99.9 \times 10^9/L$	Less than $\pm 0.5 \times 10^9/L$ Less than $\pm 5.0\%$	≥ 0.990
RBC	$0.30 \times 10^{12}/L \sim 1.00 \times 10^{12}/L$ $1.01 \times 10^{12}/L \sim 7.00 \times 10^{12}/L$	Less than $\pm 0.05 \times 10^{12}/L$ Less than $\pm 5.0\%$	≥ 0.990
HGB	20g/L ~ 70g/L 71g/L ~ 240g/L	Less than $\pm 2/L$ Less than $\pm 3\%$	≥ 0.990
PLT	$20 \times 10^9/L \sim 100 \times 10^9/L$ $101 \times 10^9/L \sim 999 \times 10^9/L$	Less than $\pm 10 \times 10^9/L$ Less than $\pm 10.0\%$	≥ 0.990

NOTICE: This brochure is for reference purpose only and subject to change without notice.

MACCURA BIOTECHNOLOGY CO., LTD.

Add: 8#, 2nd Anhe Road, Hi-tech Zone, 611731 Chengdu
Tel: +86 28 8173 1888
E-mail: maccura@maccura.com www.maccura.com



WEBSITE



WECHAT

2021-05-06-F560-2101

TÜRKİYE YETKİLİ DİSTRİBÜTÖRÜ

İZ tıbbi cihazlar İZ TIBBİ CİHAZLAR
MED. SAN. VE TİC. LTD. ŞTİ.
Egemenlik Mah. 6087 Sk. No:3 Kat:1 D:11
Tel: 0 (232) 502 58 58 • Bornova/İZMİR
info@iztibbi.com • www.iztibbi.com



İZ TIBBİ